



S7: Controls for the Duolink proximity ligation assay. Representative images from the positive and negative controls of the proximity ligation assay in figure 7. Assay performed on granulosa cells isolated from PND4 mouse ovary, red dot indicates proteins are interacting (highlighted by white arrows). (A) Positive control used SMAD3 and FOXL2 antibodies to demonstrate a confirmed interaction. (B) negative control used protein of interest (FRZB) and male-specific protein YBX2 to demonstrate no interactions. Images taken at 40x magnification. Scale bar represents 50 μm .